AKTIVLINE® BASIS

Komplet dynamisk løsning
Dynamic seating system

Dynamisk siddestilling
Dynamic seating

interco
Seat up your Life!

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Standardized dynamic seat and back including wheelchair

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* © interco GmbH, Eitorf 2011, alle rettigheder forbeholdes

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Adapter for seating system

The AKTIVLINE® basic model is a standardised care system with all proven dynamic components. It is available in 3 sizes and grows with the patient and can also be adapted to the individual patient. The width and depth of the seat and the height of the back are all designed to have 5 cm adjustment intervals, so that they can be simply and rapidly adjusted to the patient's measurements. The novel abduction wedge foot bank system (AFS)* is included in the basic configuration. The integrated abduction wedge, including the dynamic foot bank system, can be folded under the seat unit at any time and ensures an easy transfer of the patient. The entire care unit can be simply removed from the undercarriage and its weight is optimised.

- ▶ Dynamic seat system in 3 sizes
- ▶ AFS system* (abduction wedge-foot bank system)
- ▶ A variety of additional standardised options
- ▶ Can be combined with Permobil electric wheelchairs C300 and C400 (in preparation)
- ▶ Small sizes can be combined with the light undercarriage MINY®

* © interco GmbH, Eitorf 2011, All rights reserved

Medicinske krav Medical requirements

- ▶ Infantil cerebral parese: spastisk tetraplegi, dyskinetisk tetraplegi, tetraplegi med dystoniske spasmer, atetose og chorea-atetose
- ▶ Cerebral bevægelsesforstyrrelse

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- ▶ Infantile cerebral paresis: spastic tetraplegia, dyskinetisk tetraplegia, tetraplegia with dystonic spasms, athetosis and chorea-athetosis
- ▶ Cerebral movement impediments.

With this group of patients, sitting as well as the ability to make targeted movements of the head and the upper extremities is very severely restricted as a result of the following patterns:

- ▶ Extension reactions
- ▶ Spasms during extension
- ▶ Pathological movement patterns in global extension
- ▶ Episthotonus
- ▶ Startle reaction
- ▶ Fright reaction
- ▶ Synkinetic extension movements of the trunk, of the upper and lower extremities
- ▶ Dystonic spasms
- ▶ Asymmetric and symmetric tonic reflex of the back of the neck
- ▶ Tonic labyrinthine reflex


Care through the provision of a conventional, rigid seat shell is frequently not possible or very difficult to achieve, due to a lack of acceptance. In particular, the high pressure forces, resulting from conventional sitting care, are impossible for the patient to bear.

In addition, conventional sitting care increases the risk of the occurrence of: hip luxation, painful hips, shortening of the muscles, secondary damage as a result of immobility, muscle contractions or reduced mobility of the joints.

The dynamic seat shell of the AKTIVLINE seating system enables automatic changing between various body postures and thus encourages breathing, the cardiovascular system, digestion, mobility of the body joints and helps prevent the development of pressure sores.

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

- ▶ Seat back unit with physiological fulcrums
- ▶ Seat width, seat depth and back height each adjustable in 5 cm intervals
- ▶ AFS system* (abduction wedge foot bank system) in continuous design with physiological pivot points, swing-away mechanism under seat
- ▶ Anatomical basic upholstery for seat and back
- ▶ Trunk guide pads of adjustable height and width
- ▶ High quality covering in Trevira CS and spacer fabric
- ▶ Aluminium undercarriage MINY® in three seat widths (SB 35, SB 40, SB 45)
- ▶ Seat tilt with gas pressure spring by about 35°
- ▶ Seating system removable and foldable
- ▶ New tilt stop system
- ▶ Seat height and seat depth adjustable
- ▶ Wheels (12", 16", 22", 24") with drum brakes
- ▶ Push bar adjustable in height and angle

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What does “dynamic seat shell” mean?

If we observe a group of children sitting during a school class or in a cinema, we realise that the term “sitting quietly” is more of a wish, which in no way corresponds to re- ality. Everybody needs some movement even when sitting, in order to improve their own posture, to relieve parts of the body subjected to an excessive period of strain and to keep body joints mobile. AKTIVLINE is the intelligent solu- tion for all children and adults with a cerebral movement impediment, who cannot control important movement mechanisms independently. Through the dynamic linking of the various parts of the seat unit, it offers correspond- ing movement of the seat shell. This mobility is adapted to the individual requirements of the user by means of pull and push suspension, with the result that a good sitting position in the interests of therapeutic care in no way interferes with pleasant sitting.

af Pasquale Incononato, fysioterapeut med fysioterapi for revalidering, Berlin
by Pasquale Incononato, physiotherapist and rehabilitation therapist, Berlin

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What advantages does AKTIVLINE offer the user?

Many people with problems in terms of targeted movement or whose posture is influenced through mainly uncontrollable extension movements, can only benefit from care with the AKTIVLINE seat shell. In the past, attempts were made to control the involuntary movements, which for various reasons had a negative effect on sitting, by means of firm fixing. However, this type of solution only made the problem worse. The muscle tension, in particular in the hips and trunk, even increased. The consequences are common knowledge: the seat shell often broke or became deformed, the child slid out of the seat shell, the use of wedges or seat pants between the thighs caused increased pressure on the muscles, pressure which is responsible for the luxation of the hips. AKTIVLINE aims to remove these undesired consequences, and even to offer the possibility of, in part, controlling the movements described, or of using them for a sensible purpose. Thanks to the pull and push suspension, the involuntary movements are permitted to a limited extent. Once the muscles have returned to a relaxed condition, the seat shell turns back to the initial position by itself. Simultaneous physiological and dynamic sitting is achieved. Additionally, the resistance to extension does not - as is the case with conventional forms of care - take maximum effect immediately, but rather not until the end of the movement.

This means that the muscles do not develop maximum force when in a shortened condition (risk of contraction). The result: the perception of one's own body posture is noticeably improved; the possibility exists of in part internalising and controlling it.

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Which children and adults will benefit most from this type of AKTIVLINE care?

Movement impediments as a result of spastic paralysis (infantile cerebral paresis) can be observed in a wide variety of forms: on the one hand, various parts of the body can be affected (hemiparesis, diparesis, tetraparesis) and, on the other hand, movements themselves can be deficient (akinetic and hypokinetic) or not targeted (diskinetic). The perception component is always of major importance for people in the organisation of their movements. Children with tetraplegia are frequently incapable of controlling their posture and quickly become "frightened" as a reaction to various forms of irritation. An alternative form of care should be tried if over-extension of head, trunk and hips is observed, if movement of the hips is more or less bearable and if almost exclusive use is made of a wheelchair.

Nevertheless, even children who are unable to sit for more than a short time as a result of hip pain, will also feel benefits. For diskinetic movement impediments, AKTIVLINE offers the only alternative to conventional care. The main advantage lies in the fact that movement components, which disturb target-oriented actions (e.g. speaking, operating a computer), are absorbed by the movement of the shell.

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I am not happy with my rigid seat shell!

Naturally, one should not change satisfactory care through a rigid seat shell under all circumstances. It is known just how long the route to good care can be in some cases. However, if you are not satisfied despite long years of experimenting, because the seat shell frequently breaks, a good sitting posture is not achieved, because pain shortens the sitting period or because involuntary movements have a negative effect on the operation of various devices, then a trial with AKTIVLINE is a recommendable step. It will very soon become apparent whether the advantages of the “dynamic care” are accepted by the user and his/her environment.

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Status as per August 2014
Subject to technical alterations
without prior notification

Certificeret ifølge
Certified as per
DIN EN ISO 9001:2008
Reg.-nr. 73 100 357-1
DIN EN ISO 13485:2007
Reg.-nr. 70 105 357-1
TÜV Hessen