AKTIVLINE Standard C €

Monterings- og betjeningsvejledning Assembly and operating guidelines

Vision

Projektion

Konzeption

Interaktion

Innovation

Rehabilitation

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Contents

1.	Introduction	3
2.	Basic Features	4
2.1.	Adapter System	4
2.2.	Restraint System	4
2.3.	Headrest	6
2.4.	Knee Support	7
2.5.	Footrest System	8
3.	Basic Features: Undercarriage	9
3.1.	Adapter System	9
3.2.	Seat Tilt Mechanism	10
3.3.	Push Bar	11
3.4.	Anti-tip Wheels	12
3.5.	Parking Break Instructions	13
3.6.	Rear Wheels	14
3.7.	Undercarriage Folding Mechanism (optional)	15
3.8.	Pressurized Shock Absorber Function	16
3.9.	Adjustment and Assembly	17
4.	Care and Service	19
5.	Technical Data	20
5.1.	Maximum Load Specifications	20
5.2.	Tire Pressure Guidelines	20
6.	Security Tips and Recommendations	21
7.	Inspection	25
8.	Warranty	26
9.	Declaration of Conformity	27

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1. Introduction

Dear Customer,

We are very pleased that you have chosen the high quality AKTIVLINE Dynamic Seating System from interCo.

AKTIVLINE is a completely new care concept for cerebral and movement impaired individuals. The dynamic attachment of the seating system to the undercarriage construction allows for the upward movement of the seating unit. This movement potential is important, for example with the typical muscle tension and extension associated with spasticity, because it allows the seating unit to compensate by opening upward with the movement. This innovative system reduces the stress and strain associated with uncontrolled movement. Additional positive effects on the patient can be observed in patients using this new dynamic care system.

AKTIVLINE is exceptionally more robust than the normal rigid variety system. Due to the built in system dynamic, the chair materials have a longer durability.

We hope that your AKTIVLINE brings joy to your everyday activities.

The following user manual is filled with tips and insights to help make your AKTIVLINE simple and safe to use.

If you have questions, comments or suggestions regarding AKTIVLINE, please do not hesitate to contact us.

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


2. Basic Features

After the necessary fittings have been completed, you will receive your AKTIVLINE fully assembled (usually from your authorized dealer). In order to ensure that AKTIVLINE is a reliable partner for you, please pay careful attention to the following points.

2.1. Adapter System

The seat shell system is mounted to the undercarriage with an adapter system. (The seat shell undercarriage for AKTIVLINE is similar to that of the interCo SIMPLY.) This is comprised of a trapezoidal adapter (Photo 1) mounted to the underside of the shell, and the trapezoidal adapter receptacle (Photo 2) which is integrally connected to the undercarriage. The adapter is adjusted with two locking pins (Photo 3) which attach to the adapter receptacle from below and lock the adapter system securely in place.

Make certain that both locking pins and the adapter system are secured tightly.

If you determine that the locking pins have become loose, please tighten them immediately with the appropriate tool in order to guarantee the secure attachment of the shell to the undercarriage.

2.2 Restraint System

AKTIVLINE is fitted with a pelvic safety belt (Photo 4) which is essential to the correct care of the user and plays an important role in the system's function. Only when the pelvic belt is correctly fastened does the patient have the ability to utilize the AKTIVLINE movement potential and benefit from the desired effects of the system.

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The pelvic belt is secured to the shell on the right side (from the patient's perspective). On the left side (from the patient's perspective) is a buckle which allows the desired or tolerated pressure to be applied to the pelvic region.

When attaching the pelvic belt, please observe the following instructions:

The patient must be sitting as deeply as possible in the seat shell when attaching the belt. The lower back should be in direct contact to the back of the seating shell. This is most effectively accomplished with the help of the seat tilt mechanism (see Section 3.2.: Seat Tilt Mechanism) whereby the undercarriage can be tilted to the extent that the patient can be comfortably directed into the desired position.

When using a pelvic brace, please take the following precautions:

The pressure pads of the pelvic brace must lie on the two upper, front iliac spines. Be careful that the pressure of the pelvic brace does not exert pressure on the stomach or other soft tissue areas!

Feed the safety belt through the buckle of the pelvic brace and pull it tight enough that the patient sits firmly in the described sitting position (with the lower back pressed against the seat back). If this is uncomfortable for the patient, the pelvic brace should be loosened.

Keep in mind that the seat shell can open more easily when the pelvis is in continuous contact with the pelvic brace. It may take time to work up to tolerating the amount of pressure necessary to fitting the belt properly. This is accomplished by beginning with less pressure and gradually increasing the pressure until the belt lies firmly and correctly.

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2.3. Headrest

After the patient is secured in the seat with the pelvic brace, adjust the headrest. The holding devise for the AKTIVLINE headrest is integral with the back of the seat shell. To adjust the height of the support, loosen the adjustment lever at the top of the head support positioning bar (Photo 5). The angle of the headrest is adjusted with the positioning cone on the headrest bar (Photo 6).

Tip: When using the adjustment lever and there is not enough room to turn the lever a full 360°, the following trick can be used: Turn the adjustment lever as far as possible in the desired direction and then pull the lever outward. This movement unlocks the adjustment lever catch allowing the lever to be turned back without changing the selected adjustment grade.

The following points are important when adjusting the headrest:

1. A standard shaped headrest should be adjusted so that it fits directly behind the patient's head.
2. As a rule, shell shaped headrests should sit lower on the head than standard shaped supports. That means that the under side of the shell shaped headrest should ideally sit at the base of the head, placed low enough to provide support for the neck.

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2.4. Knee Support

In conjunction with the pelvic belt, the knee support (Photo 7) is another important element for the correct functioning of the AKTIVLINE. With the help of the knee support, the patient has the possibility to utilize the dynamic footrest system. The footrest dynamic is supported by the use of the foot spats, which will be discussed in the following section.

The following points are important when adjusting the knee support:

1. The knee support should be adjusted to sit directly beneath the knee cap on the patella tendon. As the knee moves up and down when the leg is bent and stretched, it is important that the knee joint is not impeded by the knee support in a way that could inhibit movement.
2. The knee support should lie on the patella tendon, but not create pressure on the area.
3. It is important that the knee support locking mechanism is firmly set in order to guarantee the correct placement of the support. Depending on how your AKTIVLINE is equipped, the knee support is either flexible with a securing strap (securely mounted to one side of the shell and attached on the opposite side with reverse rolls, padding and Velcro), or affixed with an adjustable panel and hinges (one side securely mounted, the other operable with a snap closure). With the adjustable panel variation, the knee support can be adjusted in height and depth to fit the patient.

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2.5. Footrest System

The AKTIVLINE foot rest system must be fitted to the patient's lower thigh. When doing this, it is important that the patient's feet rest comfortably on the footrest.

To adjust the footrest system to the length of the patient's lower thigh, remove the screw on the adjustment plate of the footrest and move the footrest system higher or lower along the middle bar. Once the correct position has been located, tightly reattach the screw at the correct adjustment level.

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3. Basic Features: Undercarriage

In order to familiarize you with the large variety of settings and adjustment variations, a description of the main design points of the AKTIVLINE undercarriage follows.

If uncertainties, problems or questions should arise, you are best advised to contact your authorized interCo dealer or representative in order to guarantee the safe and efficient operation of your AKTIVLINE.

3.1. Adapter System

The specifics of the adapter system are described in chapter 2.1. It is important to repeat that a safe and tight connection of the seat shell system to the undercarriage is critical.

Tip: The trapezoidal adapter receptacle is equipped with two locking pins. It is critical that the locking pins be correctly placed and secured. Before setting the seat shell on the undercarriage, both locking pins should be placed so that both setting pins project through the adapter plate. When the adapter, which is mounted under the seat shell, is now inserted into the adapter plate, the locking pins will snap into place in the adapter receptacles with an audible click sound. Test the fit of the seat to the shell system to ensure a secure and tight fit. Do this by rattling the shell and making sure that the locking pins do not open or loosen. If the shell moves or the pins open, check and readjust the adapter system. The attachment of the seat shell to the adapter should be carried out by a trained interCo representative or your authorized dealer.

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To remove the seat shell, release both locking pins by pulling in an outward direction while simultaneously turning the pins 90°. This will put the locking pins in the 'open' position. Remove the security pin which connects the back unit of the seat shell to the undercarriage. You will now be able to move the seat shell system forward and detach it from the adapter. Before reinstalling the shell, remember to turn the locking pins 90° again in order to reset them in the 'closed' position. Be careful that the security pin connecting the back unit with the undercarriage is again moved into place and that it is correctly set and adjusted.

3.2. Seat Tilt Mechanism

The AKTIVLINE undercarriage is outfitted with a mechanism allowing for the individual adjustment of the seat shell angle. Depending on the model, a continuously variable 30° angle adjustment is possible. Correcting or changing the seat angle is done with the help of two gas-pressure springs.

The following points are important when using the seat angle adjuster:

1. Press the two small release levers (Photo 8) mounted on the push bar. Allow the shell system to recline to the desired angle; then release the levers. The seat shell is now adjusted to the desired angle. If the shell does not recline to the required angle, you can assist it by applying light pressure to the push bar. To secure the position, it is important to observe that the release levers have been fully returned to their original position.

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2. By pressing the release levers, the seat angle adjuster will be released. The shell system can then be returned to its original position, assisted by the two gas-pressure springs.

3.3. Push Bar

The AKTIVLINE push bar offers two adjustment possibilities to make pushing the AKTIVLINE more comfortable:

1. First, the height of the bar is adjusted by pushing the tension lever on the back tube of the AKTIVLINE to open the mechanism. Release the lever by applying light pressure with the thumb, whereby freeing the adjustment spring (Photo 10) and allowing the push bar to be positioned at the height most comfortable for you. The push bar position will be locked into place by closing the adjustment spring.
2. In addition to varying the height of the push bar, the angle of the bar can also be individually adjusted (Photo 11). To adjust the angle, release the two clamp levers at the structural tubing joints until they are open enough to allow for the easy positioning of the bar angle. When closing the clamp lever, do not forget to ensure that the structural tubing elements are nested firmly in each other allowing a firm and secure hold. If the connection is not sufficient, it is possible that the tubing elements will slide apart leading to unnecessary structural wear and tear.

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3.4. Anti-tip Wheels

To guarantee AKTIVLINE's stability, the undercarriage is equipped with two anti-tip wheels. These are mounted with a swivel version to make curbs, stairs and sharp turns easier to negotiate.

Tip: Pay particular attention that the anti-tip wheels and the two small wheels are all pointed in the direction of the road surface when a patient is sitting in the seat!

To swivel the anti-tip wheels proceed as follows:

1. Employ the parking brake.
2. Stand behind the AKTIVLINE.
3. Pull one anti-tip bracket gently out of its tubular mounting (Photo 12).
4. Pull the anti-tip bracket outward and, in the extended position, turn it 90° so that the small wheel faces upward (Photo 13).
5. Let the anti-tip bracket spring gently back into the tubular mounting (Photo 13).
6. Repeat this procedure with the second wheel.
7. Release the parking brake.
8. Swivel back to the original position in the same way.

Tip: Pay particular attention that both anti-tip wheels always face in the same direction. The effectiveness of the anti-tip system depends on correct wheel placement.

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3.5. Parking Brake Instructions

Parking brake operating instructions:

1. Press the two brake handles on the push bar. Be careful that the brake handle does not return to its starting position (the adjustment is audible); only then is the parking brake activated.
2. Check the brakes by trying to push the AKTIVLINE forward.

Attention: If the brake does not adequately hold, bring your AKTIVLINE to an authorized dealer service center for adjustment, or repair them yourself as described in chapter 3.5.2.!

The parking brake can also be used to slow the AKTIVLINE in order to negotiate steep or hilly terrain more easily. In such cases, pull the parking brake as required, but only as far as the first or second adjustment position.

3.5.1 Parking Brake Release

Instructions for releasing the parking brake:

1. Unlock the brake from its position by releasing the adjustment lever and letting the brake lever slide back to its original position. A light pull of the adjustment lever may be necessary to release the brake lever, before the brake can be released.

3.5.2 Adjusting the Parking Brake

Tool: Combination Wrench, Size 8

Requirement: weakness in the parking break; the AKTIVLINE can still be moved despite fully applied brake.

Instructions for adjusting the parking brake:

1. Make certain that the parking brake is not applied; if necessary, release the parking brake.
2. The parking brake Bowden cables are threaded through the adjustment screw in front of the brake anchor plate, allowing adjustment of the braking power.

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3. By turning the screw counter clockwise, the brake cable is pulled tighter thereby activating the brake. The adjustment screw is then secured with a nut.
4. Remove the security nut before adjusting the brake by turning it to the left (depending on the AKTIVLINE model, a size 6 combination wrench may be required).
5. Turn the adjustment screw counter clockwise until it becomes more difficult to turn the wheel. To check this, lift the AKTIVLINE on one side and turn the effected wheel with the free hand.
6. Now turn the adjustment screw clockwise until the wheel in question again turns freely.
7. Turn the counter nut (located on the thread of the adjusting nut) clockwise until completely tightened.
8. Repeat this procedure as needed for the other wheel.

3.6. Rear Wheels

The rear wheels of the AKTIVLINE undercarriage are outfitted with quick release axles. This makes it easy to remove the wheels in case of a defective tire, or a tube needing to be replaced.

Instructions for removing or attaching the rear wheels:

1. To open the quick release wheel lever, push and hold the release button in the middle to the wheel hub while removing the wheel from the wheel receiving rod.
2. When mounting the wheel onto the receiving rod, continuously press the release button until the wheel is in place.
3. It is important to note that whether mounting or removing a wheel, the brake drums must be released in order for the wheel to either be removed or mounted.

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Tip: It is advisable to swing the anti-tip bracket downward before removing the rear wheels. This will prevent the AKTIVLINE frame from coming in direct contact with the ground, and avoid resulting damage to the carriage finish.

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3.7. Undercarriage Folding Mechanism (optional)

If you have decided on an AKTIVLINE model with a removable seat unit, it is equipped with a folding mechanism to make transport easier by minimizing the undercarriage transportation dimensions.

When folding AKTIVLINE for transport, proceed as follows:

Remove the seat shell unit from the undercarriage (to remove the shell from the adapter plate, the security stop must also be removed, since this connects the seating shell back unit with the undercarriage (Photo 15).


1. The lower ends of the back tubes are connected with a bar. The bar is in turn connected to the release with two springs which comprise the folding mechanism. Move the bar by lightly pressing down with your foot, and fold the AKTIVLINE back unit downward toward the undercarriage.
2. Match the length of the handle so that its upper side, with the help of the adjustable joints, can be folded down in front of the undercarriage front wheels.
3. Remove the rear wheels of the undercarriage from their quick release axles.

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3.8. Pressurized Shock Absorber Function

The pressurized shock absorber is designed to regulate the dynamic movement of the seat. Two adjustment screws (Photo 16) are attached to the pressurized shock absorber and enable the tension of the absorber to be set or adjusted.

- By opening both adjustment screws (turning to the left), the seat system will move and expand with the patient's uncontrolled muscle movement, tension or spasm, and then return to its original position once the movement has ended.
- By closing the upper adjustment screw (turning to the left), the dynamic seat function is deactivated and the patient will sit passively in an inflexible seat system.
- By closing the lower adjustment screw (turning to the right), the stretching of the seat shell is possible, but the seat will not automatically return to its original position. By slowly opening the lower adjustment screw the seat shell will slide back into its original position.
- In order to familiarize yourself with the various positioning possibilities, experimenting with adjusting the pressurized shock absorber will help you understand the full range of movement and possible variations.

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3.9. Adjustment and Assembly

Every customer has individual wishes and needs. These help to determine which wheel chair variation is right for them, which accessories are necessary and how the adjustments to the wheel chair and its handling should be calibrated to best suit their requirements. Working together with your physical therapist or authorized interCo dealer will assist you in finding the correct adjustments and positioning for your wheel chair.

All technical adjustments should be made only by an authorized interCo dealer or the manufacturer!

Adjusting the Seat Height and Angle:

The higher your drive wheel is attached to the frame, the steeper the seat will incline backward. This has two effects: first, the wheel chair will tilt slightly backwards, and secondly, one will sit deeper and therefore more firmly in the chair.

In combination with the inclination adjustment of the swivel wheels, the seat height can also be adjusted to optimally meet your individual requirements.

Important: By altering the drive position, the angle of the swivel wheel axle to the ground is also changed. A 90° angle should always be maintained. The knee brake must also be readjusted accordingly.

To change the installation height of the driving wheels, proceed as follows:

1. Secure the wheel chair on an even surface to avoid unintentional tipping.
2. There should not be a patient in the wheel chair.
3. Release the wheel brake.

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
Sørg for, at begge sider er indstillet ens!

4. Press the quick release axle button (on the front face of the axle) and pull the axle out.
5. Remove the drum brake security screw.
6. Adjust the axle height to the desired position and turn the security screw until it is tight.
7. Remount the wheel and the quick release axle.
8. Check that all screws and bolts have been correctly replaced and tightened, and that the brakes are functioning correctly!

Once the driving wheels have been mounted in the most comfortable position for you, the steering head inclination must be readjusted.

To guarantee perfect steering, the front wheel steering axle (dashed line) should always be perpendicular to the underground. If this is no longer the case due to alteration of the driving wheels, the front wheels can be repositioned by turning the external nuts on the front wheel receptacle (arrows). The use of a t-square (back square) is recommended to retighten the screws.

It is important that both sides are positioned identically!

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4. Care and Service

The AKTIVLINE is covered with a robust plastic coating. Most parts can be cleaned with soapy water and a soft sponge. Avoid aggressive or abrasive scrubs that can damage the plastic coating finish. All disinfection fluids standard to this area of use can be employed. The fabrics can be treated with a mild detergent.

The AKTIVLINE is mostly maintenance free. The only point requiring your regular attention is the tire pressure. The tires are therefore produced with standard air valves, so that the tire pressure can be easily tested and adjusted at your local gas station. The maximum tire pressure specification listed in the 'Technical Data' section should be adhered to and in no case exceeded.

Tip: In addition to the tire pressure, you should regularly check that all connections, screws, bolts and nuts are properly and securely fastened and that they have not become loose. Due to the demanding nature of the required care, it is not out of the question that some of the screws and fasteners will require tightening from time to time. In this case, as when brake adjustment is necessary, turn to your authorized interCo dealer. They will offer you the assistance you require or put you in direct contact with the manufacturer.

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5. Technical Data

5.1. Maximum Load Specifications

Important: The maximal load for the AKTIVLINE (working load) is 70 kg.

5.2. Admissible Tire Pressure

Tire Type	max. Tire Pressure [bar]
Front Swivel Wheels 200 x 50 mm	2,5
Rear Wheels 12"	2,5
Rear Wheels 16"	2,5
Rear Wheels 22"	4,0
Rear Wheels 24"	4,0

Important: Be careful to select the appropriate tire pressure. This decision is a compromise between the rolling comfort as experienced by the patient (lower tire pressure) and the pushing comfort of the accompanying person (higher tire pressure). Finding the right tire pressure to satisfy both needs may require a bit of experimenting. The factory suggestion is 2.5 bar.

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6. Safety Instructions

Please pay extra attention to the following safety instructions:

Conditions for transporting the AKTIVLINE-Dynamic Seating System in motor vehicles designed for handicapped persons include the use of an approved and registered wheel chair restraint system by the transporting company.

1. The maximum AKTIVLINE load is 70 kg.
2. Do not use the footrest for support when getting into or out of the wheel chair.
3. Test the effect of shifting weight (and the center of gravity) on the wheel chair; for example when driving up or down hill, on angled surfaces or when negotiating obstacles and obstructions with safety assistance from a helper. The use of an anti-tip system is strongly recommended for an inexperienced wheel chair user.
4. The seating system should be parked exclusively on horizontal, even surfaces. If parking the system in an uphill situation becomes unavoidable, be careful to bring the seat to an upright position. There is a danger of tipping backwards when the patient is in a reclined position.
5. Be careful that the seat inclination is only activated when the anti-tip lever is open.
6. Exercise caution when reaching for objects (lying in front of, behind or to the side of the wheel chair), that the patient does not lean too far out of the seat. This weight shift can create a tipping or overturning danger.
7. The anti-tip system is a device designed to prevent an unintentional backward fall. The anti-tip system should never be utilized to control the transportation function of the chair; for example, in an attempt to transport a patient when the driving wheels have been removed.
8. Stair should only be climbed with the assistance of an accompanying person. If ramps or elevators are available, they should be utilized. If ramps or elevators are not available, two people are needed to carry the

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wheel chair and patient up the stairs. If only one person is available, the anti-tip system must be activated before ascending so that the chair can not be set down on the stairs, which would otherwise lead to a dangerous fall. Afterwards the anti-tip system must be returned to the correct position. When using a loading lift with the anti-tip system, it is important to place the anti-tip device within the standard area of the loading lift.

9. The assisting person should be careful to touch the wheel chair only on the fixed, structural elements of the chair (not, for example, on the foot notch).
10. Be certain that the clamp lever of the height-adjustable push bar is closed tightly.
11. Secure the seat system undercarriage on uneven surfaces or when getting in or out of the chair (for example, into an automobile) with the brakes.
12. The knee brake, as well as the general driving handling of the AKTIVLINE, is dependent on the tire pressure. With correctly inflated drive wheels with the same air pressure on both sides, the chair will maneuver more easily and efficiently. Before driving, make certain that the air pressure is correct for your tires. The correct tire pressure is printed on the wheel cover, or should be taken from the Tire Pressure Guidelines in Section 14. The knee brake is only to be used with sufficient air pressure and the correct adjustment (max. 10mm gap, subject to technical changes).
13. It is important that the Bowden cables, which operate the drum brakes, and the gas-pressure spring do not protrude into the area of the spokes.
14. Make certain that your tires have sufficient profile.
15. We strongly warn that the person sitting in the AKTIVLINE seating system must be secured with the safety belts at all times.
16. With certain combinations or adjustment variations it is possible for the swivel wheels to collide with the footrest. If you notice that the swivel wheels are

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prohibited in their turning ability, the adjustment position should be corrected.

17. Pay close attention to traffic rules when using the AKTIVLINE on public streets and roads.
18. Wear light colored clothing in order to be more visible at night. We suggest the addition of a light to your wheel chair.
19. For upper leg amputees, the driving wheels should be installed slightly more toward the back. Activating the anti-tip system is recommended.
20. To avoid hand injuries, do not reach between the driving wheel and the knee break when rolling the chair.
21. Be careful to avoid finger injuries when using, adjusting or assembling the AKTIVLINE.
22. Retighten screws, nuts and bolts after making mechanical adjustments.
23. Take special care that the seating system is tightly affixed to the undercarriage.
24. Contact your authorized dealer or interCo representative when you have questions or require assistance with your AKTIVLINE.
25. Check the mechanical and technical connections on your AKTIVLINE regularly. Consult your authorized dealer if you are uncertain about any mechanical or technical details.
26. Remember that the AKTIVLINE is a custom product intended for the use of one individual. The individual's name is printed on the AKTIVLINE product label. The use of the AKTIVLINE for individuals other than the person for which it is designated is not permitted.
27. Attachment of additional parts or alterations to the AKTIVLINE which have not been carried out by the manufacturer is forbidden and may be legally prosecuted should it become known. In such cases, all liability claims against interCo are nullified.
28. The AKTIVLINE is intended only for use as directed. Avoid actions like driving against obstacles (such as curbs and steps) without setting the brake or jumping from elevated terrains or surfaces.

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29. Exercise caution with fire (for example, a lit cigarette)! The seat and backrest upholstery could catch fire.
30. Never leave the patient unattended; not even when the patient is secured by a seat belt and all brakes are locked.
31. Do not allow children to play with the wheel chair or its parts.
32. When the wheel chair is in motion, or when the seat and back support element are in use, other children should not play next to, under, or in front of the chair.
33. In order to eliminate injury to others, do not leave the AKTIVLINE and its passenger in the care of untrained individuals.

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7. Inspection

To ensure the optimal function of the AKTIVLINE over a long period of time, and to avoid costly repairs, we recommend a cycle of annual inspections. The cost of the annual inspection is covered by your health care organization and should only be carried out by the manufacturer.

Ask your authorized dealer about our service and maintenance contract.

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8. Warranty

In addition to our General Terms and Conditions, your AKTIVLINE product comes with a two year warranty from the date of purchase.

Defective or faulty parts will be repaired or replaced within the warranty period.

The warranty excludes modifications made during the course of development.

We accept no liability for damage caused to our product in conjunction with products of any kind manufactured by other companies, unless these have been expressly approved by interCo. Under certain circumstances, such combinations can instigate considerable damage.

Also excluded from warranty are faults due to natural wear and tear, excessive stress or loading, improper use and excessive force.

The warranty becomes invalid when original interCo parts are not used in repair work, and when the repairs or changes are not executed by interCo.

The attachment of accessories, parts or changes to the AKTIVLINE-Dynamic Seating System by anyone other than interCo will not be honored or covered by warranty, and could be liable to legal action if made known. In such cases, all liability claims against interCo are nullified.

Continuous use of your AKTIVLINE can lead to loose screws, bolts, nuts and spokes. Tighten these items regularly, or arrange a regular check up with your authorized dealer or interCo.

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*frivillige oplysninger

Guarantee postcard

If you send us this card, we can assure you an extended guarantee for your AKTIVLINE for a total of 2 years. Please fill out this guarantee postcard carefully and completely and return it to us – thank you.

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© interCo GmbH, Eitorf, 2011