# AKTIVLINE®

Den Dynamiske siddeløsning  
Dynamic seating system

Dynamisk siddestilling  
Dynamic seating

**interco**  
Seat up your Life!

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The patented, integral seating concept AKTIVLINE consists of a seat and back unit with physiological pivotal points, connected by means of sliding rails on an adapted wheelchair. AKTIVLINE offers patients physiological and dynamic seating with simultaneously low strain on body and material. AKTIVLINE contributes to both an improved sense of well-being on the part of patients as well as to their integration into society. AKTIVLINE is used wherever the previous form of care for children, young people and adults with cerebral movement impediments is unsatisfactory and problematical.

In particular children, young people and adults with infantile cerebral palsy, combined with spasms in extension or with tetraplegia and dystonic spasms, can benefit from AKTIVLINE. With patients within this framework of indications, it is frequently possible to observe that no good sitting position is achieved, resulting pain reduces the sitting time or involuntary movements render for instance eating and swallowing or the operation of certain devices impossible. Frequently, the seat shell or the wheelchair also breaks, because the forces occurring cannot be absorbed.

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### Physiological, dynamic sitting

The assumption of an orthopedically correct sitting position in line with the patient's capabilities is made possible through the individual adaptation of the seat and back units. When relaxed, the sitting angle as well as the foot and knee angles are at 90°. The anatomically-individual upholstery of the seat shell takes account of the physiology of the patient. The entire movement sequence of the shell is controlled via physiological pivotal points in the hips and knees. Seat and back units are mounted on sliding rails and enable the opening-up of the hip angle in accordance with the movement sequences of the patient. Knee-angle and footrest unit simultaneously permit the dynamic extension tendency of the patient. A readjusting spring system in both the seat as well as the leg system automatically returns the patient to the physiological starting position.

By allowing the natural urge for and sequence of movement, it is possible to keep the moving parts of the body mobile. Uncontrolled extension movements of the patient are permitted in a controlled manner. The position of head and shoulders in space will change minimally with maximal stretching of the body and chair. The orientation will be preserved. Activities like watching, drinking or operating a joystick aren't disturbed. Vestibular and visual perceptions are stable.

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### Reduced burdens

The dynamic bearing arrangement of the seat shell on the wheelchair absorbs emerging forces. In contrast to conventional, rigid seat shells with a rigid sitting position, AKTIVLINE alters the seat and leg-angle position in the event of cramp situations arising and automatically adjusts to the body extension of the patient. Emerging forces are not transmitted directly to body and material, but rather are absorbed by the admissible individual movement radius. The detailed, stable construction of seating system and undercarriage is simple to operate and can be adapted to the patient and accompanying person in many ways. Due to the comprehensive technology, the overall weight of the complete wheelchair unit is not inconsiderable; nevertheless, this can be reduced individually through the omission of various individual components. For carers and parents, however, the use of AKTIVLINE offers one decisive advantage: there is no need for permanent correction of the patient's sitting position.

### Improved sense of well-being

AKTIVLINE enables the patient to perform targeted actions such as speaking, swallowing or even operating a computer. Involuntary movements, caused for example by joy or fright as well as other disturbances of natural movement sequences are absorbed through the mobility of the seat shell, thus enabling concentration on the desired, sensible movement. Firm fixing and rigid sitting lead to increased

Individuelle behov kræver individuelle løsninger  
Individual requirements necessitate individual solutions

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muscle tenseness and thus to discomfort. The sitting period is reduced. Movements help prevent tenseness increasing yet further and thus enable longer sitting periods.

The term "sitting quietly" does not correspond to reality. Every person needs some movement when sitting, in order to improve their own body posture, to relieve parts of the body subjected to strain for too long and to keep joints mobile. AKTIVLINE enables movement when seated with simultaneous orthopedically-correct posture and contributes towards a longer sitting period. Pain as a result of tension can be reduced. Physiological and dynamic sitting conforms to modern-day medical-therapeutic requirements and is supported through the use of the integral AKTIVLINE concept. The patented construction principle of AKTIVLINE contributes to considerably reduced strain on both the patient and the material. This results in a considerable lowering of costs as well as in a fundamental improvement in the patient's sense of well-being. There are no more obstacles to social integration.

We recommend that you select the suitable technical aid for your child together with your specialist doctor, therapist and with a specialist medical equipment dealer of your choice. We at interco shall be pleased to provide you with assistance on site. Questions? Simply give us a call – we shall be glad to help.

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- ▶ Infantile cerebral paresis: spastic tetraplegia, dyskinetic tetraplegia, tetraplegia with dystonic spasms, athetosis and chorea-athetosis
- ▶ Cerebral movement impediments.

With this group of patients, sitting as well as the ability to make targeted movements of the head and the upper extremities is very severely restricted as a result of the following patterns:

- ▶ Extension reactions
- ▶ Spasms during extension
- ▶ Pathological movement patterns in global extension
- ▶ Episthotonus
- ▶ Startle reaction
- ▶ Fright reaction
- ▶ Synkinetic extension movements of the trunk, of the upper and lower extremities
- ▶ Dystonic spasms
- ▶ Asymmetric and symmetric tonic reflex of the back of the neck
- ▶ Tonic labyrinthine reflex

Care through the provision of a conventional, rigid seat shell is frequently not possible or very difficult to achieve, due to a lack of acceptance. In particular, the high pressure forces, resulting from conventional sitting care, are impossible for the patient to bear.

In addition, conventional sitting care increases the risk of the occurrence of: hip luxation, painful hips, shortening of the muscles, secondary damage as a result of immobility, muscle contractions or reduced mobility of the joints.

The dynamic seat shell of the AKTIVLINE seating system enables automatic changing between various body postures and thus encourages breathing, the cardiovascular system, digestion, mobility of the body joints and helps prevent the development of pressure sores.

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An important utensil: the concept form.

#### Basic features

The AKTIVLINE caresystem contains an aluminium seat and backrest connected by physiological pivot points and an aluminium undercarriage SIMPLY® or SIMPLY LIGHT®. Seat and backrest have a good anatomically molded groundpolster including side support for the upper body and the hips and upper legs. Different options will be chosen according to the individual needs. AKTIVLINE seat and back unit in the above basic version is for assembly on an electric wheelchair provided by the customer or an indoor undercarriage ROOMY NE incl. electric motor.

#### Product variants

The AKTIVLINE care system is available in standardized and individual version. The individual version is possible in aluminium or carbon seat and backrest.

Your personal, local specialist advisor will assist you in the planning of an individual solution, taking account of your personal requirements. We at interco guarantee the availability of replacement parts or readjustment over a period of years

#### Accessories/special equipment

A comprehensive range of accessories is available and is chosen carefully in accordance with the requirements of the patients and accompanying persons.

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### What does “dynamic seat shell” mean?

If we observe a group of children sitting during a school class or in a cinema, we realise that the term “sitting quietly” is more of a wish, which in no way corresponds to reality. Everybody needs some movement even when sitting, in order to improve their own posture, to relieve parts of the body subjected to an excessive period of strain and to keep body joints mobile. AKTIVLINE is the intelligent solution for all children and adults with a cerebral movement impediment, who cannot control important movement mechanisms independently. Through the dynamic linking of the various parts of the seat unit, it offers corresponding movement of the seat shell. This mobility is adapted to the individual requirements of the user by means of pull and push suspension, with the result that a good sitting position in the interests of therapeutic care in no way interferes with pleasant sitting.

af Pasquale Inconorato, fysioterapeut med fysioterapi for revalidering, Berlin  
by Pasquale Inconorato, physiotherapist and rehabilitation therapist, Berlin

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## What advantages does AKTIVLINE offer the user?

Many people with problems in terms of targeted movement or whose posture is influenced through mainly uncontrollable extension movements, can only benefit from care with the AKTIVLINE seat shell. In the past, attempts were made to control the involuntary movements, which for various reasons had a negative effect on sitting, by means of firm fixing. However, this type of solution only made the problem worse. The muscle tension, in particular in the hips and trunk, even increased. The consequences are common knowledge: the seat shell often broke or became deformed, the child slid out of the seat shell, the use of wedges or seat pants between the thighs caused increased pressure on the muscles, pressure which is responsible for the luxation of the hips. AKTIVLINE aims to remove these undesired consequences, and even to offer the possibility of, in part, controlling the movements described, or of using them for a sensible purpose. Thanks to the pull and push suspension, the involuntary movements are permitted to a limited extent. Once the muscles have returned to a relaxed condition, the seat shell turns back to the initial position by itself. Simultaneous physiological and dynamic sitting is achieved. Additionally, the resistance to extension does not - as is the case with conventional forms of care - take maximum effect immediately, but rather not until the end of the movement.

This means that the muscles do not develop maximum force when in a shortened condition (risk of contraction). The result: the perception of one's own body posture is noticeably improved; the possibility exists of in part internalising and controlling it.

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### Which children and adults will benefit most from this type of AKTIVLINE care?

Movement impediments as a result of spastic paralysis (infantile cerebral palsy) can be observed in a wide variety of forms: on the one hand, various parts of the body can be affected (hemiparesis, diparesis, tetraparesis) and, on the other hand, movements themselves can be deficient (akinetic and hypokinetic) or not targeted (diskinetic). The perception component is always of major importance for people in the organisation of their movements. Children with tetraplegia are frequently incapable of controlling their posture and quickly become "frightened" as a reaction to various forms of irritation. An alternative form of care should be tried if over-extension of head, trunk and hips is observed, if movement of the hips is more or less bearable and if almost exclusive use is made of a wheelchair.

Nevertheless, even children who are unable to sit for more than a short time as a result of hip pain, will also feel benefits. For diskinetic movement impediments, AKTIVLINE offers the only alternative to conventional care. The main advantage lies in the fact that movement components, which disturb target-oriented actions (e.g. speaking, operating a computer), are absorbed by the movement of the shell.

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### **I am not happy with my rigid seat shell!**

Naturally, one should not change satisfactory care through a rigid seat shell under all circumstances. It is known just how long the route to good care can be in some cases. However, if you are not satisfied despite long years of experimenting, because the seat shell frequently breaks, a good sitting posture is not achieved, because pain shortens the sitting period or because involuntary movements have a negative effect on the operation of various devices, then a trial with AKTIVLINE is a recommendable step. It will very soon become apparent whether the advantages of the “dynamic care” are accepted by the user and his/her environment.

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Status as per August 2014  
Subject to technical alterations  
without prior notification

Certificeret ifølge  
Certified as per  
**DIN EN ISO 9001:2008**  
Reg.-nr. 73 100 357-1  
**DIN EN ISO 13485:2007**  
Reg.-nr. 70 105 357-1  
TÜV Hessen

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© interco GmbH, Eitorf, 2014