



MOTIONCOMPOSITES

*The lightest
HD wheelchair
in the world!*

350 lbs maximum weight capacity - up to 22" wide

HELIO C2 HD

By using a new extra large carbon fiber weave to reinforce the Helio C2 HD crossbrace, we were able to increase the maximum weight capacity significantly while adding minimal weight to the wheelchair.

The new Helio C2 HD can handle up to 350 lbs with only a 176 gram weight increase over the standard Helio C2*.

**Heavy duty
without the HEAVY.**



**2014 AWARD
HARDING
INNOVATIVE
PRODUCT**



Canadian Seating &
Mobility Conference

The strong yet lightweight design of the Helio C2 HD earned the Harding Innovative Product Award at the 2014 Canadian Seating and Mobility Conference.

Widths 18, 19 and 20 inches available with 350 lbs crossbrace.
Widths of 21 and 22 inches only available with 350 lbs crossbrace.

*20"x20" Helio C2 HD compared to 20"x20" Helio C2.

www.motioncomposites.com